

Singfiller®

sterile injection
safety and natural
strong resistance

CE
1023




Cross-linked sodium Hyaluronate for Plastic Surgery

COMPANY PROFILE

Hangzhou Singclean Medical Products Co., Ltd established in March 2003 with registered capital of USD13,000,000 is a national high-tech enterprise integrating the R&D, production and sales of medical biomaterials. With well-qualified staff, 5S Standard and first-class facilities, the company ensures the continuous production of medical devices products of Class III that conforms to the laws and regulations of the State Food and Drug Administration (SFDA) and EU MDD.

Sticking to the quality policy of “Quality First, Management Priority, Building a Harmonious Enterprise, Creating the Singclean Brand”, the company passed the ISO9001 and ISO13485 quality management system in 2007, and obtained the Innovation Fund Certificate by Ministry of Science and Technology with the Microporous Polysaccharide Hemospheres project in 2008. The company obtained the certification of China high-tech enterprise in 2010. The medical sodium hyaluronate gels passed the CE certificate in 2011, and the cross-linked sodium hyaluronate gel for plastic surgery brand— **Singfiller**® passed the CE certificate in 2017.



WHAT IS SINGFILLER?

Singfiller is BDDE-cross-linked hyaluronic acid filler which is colorless, transparent, viscous, sterile non-pyrogenic injection. This third generation product is reticulated and obtained via bio-fermentation technology.


Singfiller is developed with cutting-edge technology, guaranteed for purity, safety and long-lasting effect.




TYPE	Fine	Derm
COMPOSITION	20mg/ml(HA)	20mg/ml(HA)
Particle size	0.10~0.15mm	0.15~0.28mm
Recommended indications	thin superficial lines ②③④⑤	moderate wrinkle ③④⑤⑥⑦⑧
Where to inject	upper part of dermis	middle part of dermis
Volume of Syringe	1mL (CE)	1mL (CE) / 2mL (CE)
Storage	0~30°C (32-86°F)	0~30°C (32-86°F)
Use	Single	Single
Duration	6 - 9 months	6 - 9 months
Needle 2.25PRTC	27G*1	27G*1
1.0PRTC	30G*2	30G*2

APPLICATION

- 1 facial beauty injection
- 2 improve skin appearance
- 3 structure, reduce wrinkle, prevent skin aging
- 4 improve labial shape
- 5 facial tissue augmentation



- ① Nasolabial folds
- ② Frown lines
- ③ Wrinkles across the forehead
- ④ Neck lines
- ⑤ Nose
- ⑥ Chin
- ⑦ Lip
- ⑧ Earlobe
- ⑨ Temple
- ⑩ Cheek




Deep	Deeper	Sub-Q
20mg/ml(HA) 0.28~0.50mm	20mg/ml(HA) 0.50~1.00mm	20mg/ml(HA) 1.25~2.00 mm
Deep facial wrinkles and folds ①⑥⑦⑧⑨⑩ deep layer of dermis	Deep facial wrinkle and folds ①⑨⑩ deep dermis or subcutaneous tissue shallow	the local depression, soft tissue filler. And breast enhance , Breast butt deep dermis or subcutaneous tissue shallow
1mL (CE) / 2mL (CE)	1mL (CE) / 2mL (CE) / 10mL(No CE)	2mL(No CE) / 20mL(No CE)
0~30°C (32-86°F)	0~30°C (32-86°F)	0~30°C (32-86°F)
Single	Single	Single
6 - 9 months	6-12months	6-12months
26G*1 27G*2	26G*1 26G*2	/ /


NASOLABIAL INJECTION



MECHANISM OF ACTION



Singfiller® shows equivalent results on Suitability, Safety & Effectivity.





Singfiller®

Cross-linked sodium Hyaluronate for Plastic Surgery



www.singclean.net
www.singfiller.org

Hangzhou Singclean Medical Products Co., Ltd
Add: No. 125 (E), 10th Street, Hangzhou Economic and
Technological Development Zone, ZheJiang, China
P.C.: 310018
Tel: +86-571-63431868 Fax: +86-571-63431886
E-mail: sales@hzxhe.com
Web: www.singclean.net www.singfiller.org